



B.BRAUN VET CARE

YOUR PARTNER IN IV THERAPY

Caresite® Needlefree Vascular Access System

- Superior microbial barrier
- Helps to prevent key-part touch contamination
- Allows to visually assess and confirm that a proper flush has been performed



CODE	DESCRIPTION	UNITS	MSRP
BB-415122-01	Caresite® Needlefree Valve	Pack 100	\$159.99
BB-470100-01	Caresite® Needlefree Single Extension	Pack 100	\$296.13

IntroVet® IV catheter The Polyurethane catheter

- Facilitates one-handed catheter advancement
- Transparent flashback chamber allows quick visualization of blood
- Permits attachment of a syringe for aspiration or other specialized procedures
- Sharp needle bevel offers a reduction in pain due to lower forces
- Creates a V-shaped, tricuspid incision versus a lancet cut for easier catheter insertion



CODE	DESCRIPTION	UNITS	MSRP
BB-4251300V	IntroVet® PUR 24G x 3/4" Yellow PUR 24 3/4" 22 ml	Pack 50	\$58.63
BB-4251318V	IntroVet® PUR 22G x 1" Blue PUR 22 1" 35 ml	Pack 50	\$58.63
BB-4251326V	IntroVet® PUR 20G x 1 1/4" Pink PUR 20 1 1/4" 60 ml	Pack 50	\$58.63
BB-4251334V	IntroVet® PUR 18G x 1 1/4" Green / white PUR 18 1 1/4" 105 ml	Pack 50	\$58.63

Intrafix® IV Administration Set

- No need to reprime the set if air enters the line
- Helps to reduce your work time
- No air in line so a new set is not required
- Helps to prevent air embolism and infusion related infections
- Helps to reduce fluid spills and needlestick injuries (Air Inlet)



CODE	DESCRIPTION	UNITS	MSRP
BB-4063000	Intrafix® SafeSet II Tube Length 70" PVC, DEHP free	Pack 100	\$195.44
BB-4063003	Intrafix® SafeSet II Tube Length 90" PVC, DEHP free	Pack 100	\$302.05

Mini-Spike® 2 Withdrawal and Injection Spike for multidose vials.

- Swabable valve automatically closes upon disconnection of syringe
- Needle-free access Compatible with Luer-Lock and Luer-Slip
- Integrated air vent making over-pressure release techniques unnecessary
- Allows multiple access for direct aspiration or fluid injection
- Two-lumen spike fitting with common vials and containers
- Designed to ensure clean, leak-proof and convenient injections and withdrawals



CODE	DESCRIPTION	UNITS	MSRP
BB-4550591	Mini-Spike® 2 Filter Swabable Needlefree Valve	Pack 50	\$177.68

Omnifix® Three-Piece Syringes For Smooth, Controlled Delivery

- Smooth plunger action and double-seal stopper
- Secure backstop to prevent accidental plunger withdrawal
- High visibility – transparent barrel with bold graduations



CODE	DESCRIPTION	UNITS	MSRP
BB-4617100V	Omnifix® Three-piece Syringes 10 ml Luer Lock	Pack 100	\$31.98
BB-4617207V	Omnifix® Three-piece Syringes 20 ml Luer Lock	Pack 100	\$69.89
BB-4617509F	Omnifix® Three-piece Syringes 50 ml Luer Lock	Pack 100	\$159.91

Extension lines T-Port Extension with injection port and side clamp

- Small bore infusion extension with side clamp and injection site
- Outer ø: 2.4 mm, inner ø: 1.3 mm, Volume 0.1 ml, Length 10 cm



CODE	DESCRIPTION	UNITS	MSRP
BB-471954	Small bore extension with T-port and spin lock	Pack 100	\$302.05

Discofix® C Stopcocks

- Allows convenient access for injections and aspirations by rotating adapter
- Lets you feel the proper setting by tactile response at every 45° turn
- Resists the hazardous effects of pharmaceuticals
- Maintains a leakproof connection even during long-term application
- Helps to prevent leakage, contamination and air embolism



CODE	DESCRIPTION	UNITS	MSRP
BB-16494C	Discofix® C 3-way Stopcock Blue	Pack 100	\$165.83

Venofix® Winged IV Needle

- Venipuncture sets for short-term injection, infusion, transfusion
- Winged IV needle with 30 cm tube and closing cone



CODE	DESCRIPTION	UNITS	MSRP
BB-4056337	Venofix® Gauge 21 Tube 12" Flow/Min 43ml Green	Pack 50	\$29.02
BB-4056310	Venofix® Gauge 19 Tube 12" Flow/Min 22ml Beige	Pack 50	\$29.02